

“Calibrated” Refractive Instruments

(PATENT PENDING)

For “Myopia, Hyperopia and Enhancements”



DAVID J. SCHANZLIN, M. D.

New “Calibrated” Laser Refractive Instruments provide the Surgeon with the added information they require during Laser Refractive Procedures.

- “Real Time, Intra-Operative Measurements means better results.”

34-690



SCHANZLIN “LASIK Irrigating Cannula”, Flattened 26g with Calibrating Scale.
(Patent Pending)

- Measures Diameters of Microkeratome Cuts for a Calibrated Overview of Laser Ablation Patterns.

34-7100



CARO “ LASIK Flap Manipulator” Curved Tip with Calibrating Scale.
(Patent Pending)

- Measures Any Flap or Corneal (Epithelial) Irregularities inside or outside of the desired Laser Ablation Patterns.

34-695



SCHANZLIN “Pickup & Utility Forceps” Cross-Action with Calibrating Scale.
(Patent Pending)

Precision, Quality, Consistency and Performance

Made in the U.S.A.

B. Graczyk, Inc. - Universal™ / Surgical Eye Care Products

61 Brookside Drive - Glendale Heights, IL 60139

Phone: (630) 871-2440 - Fax: (630) 871-7620

Email: graczykinc@aol.com Web Site: www.keratomeblades.com




"CALIBRATED" REFRACTIVE INSTRUMENTS

(PATENT PENDING)

The New "Calibrated" Laser Refractive Instrument designs provide the surgeon with the added information they require during all Laser Refractive Procedures.

"Myopia, Hyperopia and Laser Enhancements"




After using "Calibrated Scale"

"Real Time, Intra-Operative" measurements using a "Calibrated Scale" placed on the anterior surface of the instruments. (Patent Pending)

- Instruments measure the diameters of the LASIK Flap.
- Measures Diameters of Microkeratome Cuts for a Calibrated Overview of Laser Ablation Patterns prior to use.
- Measures any LASIK Flap or Corneal (Epithelial) Irregularities inside or outside of the desired Laser Ablation Pattern in "Real-Time."


34-690

Schanzlin LASIK Irrigating Cannulas & Flap Manipulator, flattened 26 gauge, angled 11mm with "Calibrated Scale." (Patent Pending)




34-7100

Caro LASIK Flap Manipulator & Elevator, curved spatula tip with "Calibrated Scale." (Patent Pending)



34-695

Schanzlin Shield Pickup / Utility Forceps, curved to flat, cross-action with "Calibrated Scale." (Patent Pending)



Instruments are Stainless Steel, Reusable, Autoclavable and made in the USA.
All Universal™ / Surgical Instruments are Guaranteed for Life.

LASIK Cal Inst / Sch