


## DUO-CELL™ PVA Bacteriostatic Products

Duo-Cell™ “FDA Approved” Patented Technology Inoculates our Hydrophilic, High Density and Lint-Free PVA formulation with “Bacteriostatic Properties”. Duo-Cell™ inhibits the growth of microorganisms; highly absorbent PVA Material traps fluids and contaminants deep inside effectively reducing the risk of Bio-Burden Cross-Contamination during “EVERY” surgical procedure.

**32-101**



**Blunt Sponge**

**32-111**



**Schanzlin  
Drain**

**32-103**



**Spear**



- Duo-Cell™ PVA is Non-Toxic and Non-Leaching
- Biocompatible
- High Density PVA Formulation
- Hydrophilic, Collects Fluids and Traps Contaminates
- Prewashed - Lint Free
- Highly Resistant to Chemicals,
- Multiple uses and Applications
- Packaged Sterile, Single Use
- Gamma Radiation Sterilized

**32-141**



**Instrument Wipes**

**32-151**



**8" x 8" x 2 mm  
Sterilization Tray or  
Mayo Stand Pads**

**32-131**



**Diamond Knife  
Cleaning Block**

***Precision, Quality, Consistency and Performance***

Made in the U.S.A.

**B. Graczyk, Inc. - Universal™ / Surgical Eye Care Products**

61 Brookside Drive - Glendale Heights, IL 60139


Phone: (630) 871-2440 - Fax: (630) 871-7620

Email: [graczykinc@aol.com](mailto:graczykinc@aol.com) Web Site: [www.keratomeblades.com](http://www.keratomeblades.com)



# DUO-CELL™ PVA SURGICAL PRODUCTS


Universal™ / Surgical Eye Care Products



32-101

**DUO-CELL™ BLUNT SPONGE (PREWASHED)**  
PVA Surgical Eye Sponge, Non-Linting, Quick Absorption with Squared Corners, Blunt Tip on a Tapered Malleable Handle.

**Sterile “Single-Use”, 5 Per Pack, 20 Packs Per Box**



32-103

**DUO-CELL™ EYE SPEAR**  
PVA Surgical Eye Spear, Non-Linting, Quick Absorption with Rounded Corners. Triangle Shape on a Tapered Malleable Handle.

**Sterile “Single-Use”, 5 Per Pack, 20 Packs Per Box**



32-111

**DUO-CELL™ SCHANZLIN “LASIK” DRAIN, Universal Design**  
Accommodates Either “Superior or Nasal LASIK Flaps. PVA Drain to Wick Excess Fluids, 11.5mm Aperture Opening


**Sterile “Single-Use”, 5 Per Pack, 20 Packs Per Box**



32-131

**DUO-CELL™ Diamond Knife Cleaning Block, Designed for**  
Cleaning and Removing Debris from the Diamond Blade.


**Sterile “Single-Use”, 1 Per Pack, 10 Packs Per Box**



32-141

**DUO-CELL™ Surgical Instrument Wipes, Designed for Lint-Free**  
Cleaning of Delicate Surgical Instruments. Absorbs Body Fluids while Cleaning the Instrument.

**Sterile “Single-Use”, 3”x 3” x 2mm, 1 per Pack, 20 Packs Per Box**



32-151

**DUO-CELL™ Multipurpose PVA Pad for Mayo Stand,**  
Sterilization Tray, Instrument Assembly or Cleaning Stations. High Density PVA Formulation, Hydrophilic, Collects Fluids and Traps Contaminates.

**Sterile “Single-Use”, 1 Per Pack, 10 Packs Per Box**

DUO-CELL™ is safe for use in all steam autoclaves.

## Ordering Information: